

Poster Reprint

AOAC 2025
Poster number W020

Method Translation for the Analysis of Vanilla Extracts Using a Compact, Single Channel GC-FID and Carrier Gas Switching Module

Jessica Kovach and Lorna De Leoz
Agilent Technologies, Wilmington DE

Introduction

Every year consumers grow more nutrient-conscious, steering their consumption habits towards ingredients labeled as natural or even organic. This trend even extends into flavors that are called out in an ingredient list on a food or beverage product.

Vanilla is a flavor variety that has been heavily impacted by changes in consumer habits as well as by supply and demand of vanilla beans used to produce high-quality vanilla extracts. Unfortunately, these market conditions, increase the likelihood of adulteration in vanilla extracts by use of aroma chemicals that are not naturally occurring. Common adulterants are coumarin (C), ethyl vanillin (EV), eugenol (E), and guaiacol (G).^{2,3} Adulteration impacts the standard of identity for vanilla extracts established by the U.S. FDA; therefore, flavor companies rely on QC labs to implement reliable and sensitive analysis techniques for screening incoming raw ingredients against expected purity standards before the ingredient can enter a flavor formula on the manufacturing floor. Screening is typically executed with gas chromatography (GC)-based methods adapted from those used in a R&D lab. However, differences in column length, diameter, and phase can make it seem challenging to harmonize GC methods across the R&D and QC functions within a company due to the different analysis objectives of the respective laboratories.

Experimental Goals-

- Demonstrate the ease of the Agilent Method Translator in harmonizing GC methods for multiple columns and gases
- Establish method precision, linearity, and resolution on standards of vanilla-like adulterants using the translated methods
- Successfully showcase the method translation using store-bought extracts to verify the integrity of the fastest GC method on samples containing matrix
- Screen for and quantify any vanilla adulterants in the samples of vanilla extracts with the goal of detecting no adulterants using the fastest GC method.

Experimental

Table 1. GC Parameters used in the Method Translator

Parameter	Value
GC System	Agilent 8850 GC with 7693A Automatic Liquid Sampler
ALS	1 μ L injection Solvent A = isooctane, 1 prewash, 1 post wash Solvent B = isooctane, 1 prewash, 1 post wash, 1 sample wash, 6 sample pumps S/SL Syringe: 10 μ L (p/n G4513-80203)
Split/Splitless Inlet	325 °C Septum purge: 3 mL/min 1) 25:1 split 2) 25:1 split 3) 50:1 split 4) 200:1 split - Inlet septa, advanced green, nonstick (p/n 5190-3158) - Low pressure drop split liner (p/n 5183-4647) - Column nut for GC capillaries (p/n 5181-8830) - Column connection - 6 mm using graphite ferrule tool (p/n G3440-80217)
Helium Conservation Module	- Output channel plumbed to S/SL EPC - AUX conservation gas channel plumbed to house hydrogen (H ₂) gas - Helium channel plumbed to house helium (He) gas
Column	1) Agilent J&W DB-1 60 m \times 320 μ m, 0.25 μ m (5' cage) (p/n 123-1062E) Graphite Ferrules 0.1 to 0.32 mm column (p/n 5080-8853) 2) Agilent J&W DB-1 30 m \times 250 μ m, 0.25 μ m (5' cage) (p/n 123-1032E) Graphite Ferrules 0.05 to 0.25 mm column (p/n 500-2114) 3) Agilent J&W DB-1 20 m \times 180 μ m, 0.18 μ m (5' cage) (p/n 123-1022E) Graphite Ferrules 0.05 to 0.25 mm column (p/n 500-2114) 4) Agilent J&W DB-1 10 m \times 100 μ m, 0.10 μ m (5' cage) (p/n 123-1012E) Graphite Ferrules 0.05 to 0.25 mm column (p/n 500-2114) Constant flow (He): 1) 1 mL/min 2) 1 mL/min (best efficiency translation) 3) 0.72 mL/min 4) 0.4 mL/min Constant flow (H ₂): 1) 1 mL/min 2) 1.25 mL/min (best efficiency translation) 3) 0.90 mL/min 4) 0.5 mL/min
Oven	He carrier parameters: 1) 40 °C (hold 0 min), 5 °C/min to 280 °C (2 min hold); run time = 50 min 2) 40 °C (hold 0 min), 12.556 °C/min to 280 °C (0.8 min hold); run time = 19.91 min 3) 40 °C (hold 0 min), 21.136 °C/min to 280 °C (0.48 min hold); run time = 11.84 min 4) 40 °C (hold 0 min), 51.36 °C/min to 280 °C (0.2 min hold); run time = 4.87 min H ₂ carrier parameters: 1) 40 °C (hold 0 min), 5 °C/min to 280 °C (2 min hold); run time = 50 min 2) 40 °C (hold 0 min), 14.756 °C/min to 280 °C (0.68 min hold); run time = 16.94 min 3) 40 °C (hold 0 min), 27.109 °C/min to 280 °C (0.37 min hold); run time = 9.22 min 4) 40 °C (hold 0 min), 69.974 °C/min to 280 °C (0.14 min hold); run time = 3.57 min * Optional: In place of post-run bakeout, add a second ramp after 280 °C to 325 °C at 60 °C/min to perform an in-run column bakeout Oven equilibration = 1 minute Post run bakeout = 325 °C for 2 minutes (* optional for column bakeout)
FID	300 °C, H ₂ = 30 mL/min, air = 400 mL/min, N ₂ = 25 mL/min Universal 0.011 inch id FID jet (p/n 5200-0176)
Data Rate	20 Hz
Injections in Sequence	Ten injections for 3 days of interday precision; helium injections followed by hydrogen injections Four injections per concentration level in linearity and working calibration; helium injections followed by hydrogen injections Six injections per vanilla extract variety; helium injections followed by hydrogen injections




Figure 1. Agilent 8850 GC equipped with a helium conservation module where the AUX channel is configured to use H₂ as a second carrier gas for sequence efficiency

Results and Discussion

Method Translation-

A multi-standard of vanillin, an expected analyte, along with the four possible vanilla adulterants was prepared at 100 ppm. The method translation achieved good chromatographic peak shape and baseline resolution (Figure. 3 and Table 2) using both He and H₂ carrier gases across four DB-1 columns of varying lengths.

Table 2. Resolution of analytes at 100 ppm during the 60 m to 10 m method translation

Compound	60 m Retention time (min)	Average Resolution H ₂ Carrier (n=10)	10 m Retention time (min)	Average Resolution H ₂ Carrier (n=10)
Decane (ISTD)	16.254	88.995	1.260	16.471
Guaiacol	18.092	25.169	1.401	18.403
Eugenol	25.574	97.835	1.962	69.526
Unknown Impurity	25.858	3.725	n.d.	n.d.
Vanillin	26.165	4.024	2.010	5.818
Coumarin	27.160	12.183	2.088	9.175
Ethyl Vanillin	27.678	6.444	2.122	3.979

* n.d. = Not detected

Interday Area Precision-

Area precision was conducted on the analytes with ten injections per day per carrier gas then were normalized to the internal standard and averaged over three days of experimentation (Table 3).

Table 3. Interday precision (%RSD) for the four J&W DB-1 columns using H₂ carrier gas

Average area precision (%RSD) (n = 10 Injections/Day)					
Column	Guaiacol	Eugenol	Vanillin	Coumarin	Ethyl Vanillin
Day 1					
60 m	2.114	1.789	1.487	1.487	1.597
30 m	0.885	0.956	1.013	0.930	0.997
20 m	0.649	0.801	0.844	0.849	0.850
10 m	1.125	0.818	0.995	0.966	1.067
Day 2					
60 m	1.429	1.537	1.424	1.595	1.594
30 m	0.707	0.676	1.227	0.732	0.970
20 m	0.467	0.560	0.646	0.509	0.621
10 m	1.119	0.920	1.096	0.908	0.984
Day 3					
60 m	1.224	1.357	2.441	1.780	2.223
30 m	0.792	1.110	1.434	1.170	1.317
20 m	0.590	0.655	0.722	0.569	0.666
10 m	0.582	0.740	0.746	0.749	0.900




Figure 2. Agilent Method Translator for the translation of a 60 m to a 10 m DB-1




Figure 3. Method translation with H₂ carrier gas, four columns and five vanilla-like analytes

Results and Discussion

Method Linearity-

Linearity is demonstrated for each of the analytes of interest from 10 to 100,000 ppm using both He and H₂ carrier gases. Demonstrating linear behavior across a wide concentration range is critical because the use rate of the adulterants can be quite large and variable; therefore, a method suitable across a large concentration range is highly critical in the flavor space. The R² ranged from 0.9997 to 1.000

Analysis of Vanilla Extracts and Flavor-

The area precision of each extract and flavor under both carrier gas conditions on a 10 m DB-1 was < 4% RSD for the detectable analyte area. It is important to assess precision of the fast analytical methods with the vanilla extracts to look for potential matrix effects that would make the translated method harder to use in the QC environment of a flavor house.

Quantification of Detectable Analytes-

To better quantify any detectable quantities of the five vanilla aroma compounds in store bought samples, a working calibration range in addition to method linearity was studied from 100 to 5,000 ppm. The tighter concentration range was chosen because it is more relevant to the expected analyte concentrations in these store-bought products. Coefficients of determination ranged from 0.9998 to 1.000.




Figure 5. Detectable analyte concentration (ppm) for three vanilla extract market products using H₂ carrier gas

Compound	R ²	
	Helium Carrier Gas	Hydrogen Carrier Gas
30-meter DB-1 (n=4)		
Guaiacol	0.9999	0.9999
Eugenol	1.0000	0.9999
Vanillin	0.9997	0.9997
Coumarin	0.9998	0.9998
Ethyl Vanillin	1.0000	0.9999
10-meter DB-1 (n=4)		
Guaiacol	1.0000	1.0000
Eugenol	0.9999	0.9999
Vanillin	0.9998	0.9998
Coumarin	0.9998	0.9998
Ethyl Vanillin	1.0000	1.0000

Table 4. Coefficient of determination for linearity ranging 10 to 100,000 ppm on the 30 m and 10 m DB-1

Conclusions

Method Translation-

- The method translator tool provided an easy conversion of a 50-minute-long method to one that is < 5 minutes resulting in a speed gain of 10-fold and 14-fold, for He and H₂ carrier gases, respectively (Figures 2 and 3)
- Resolution of analytes were >4 on the 60 m and >3.5 on the 10 m DB-1 columns (Table 2)
- Precision was <1.85% RSD and <2.5% RSD for He and H₂ carrier, respectively (Table 3)
- Linearity was demonstrated from 10 ppm to 100,000 ppm for He and H₂ carrier with R² ranging from 0.9997 to 1.0000

Carrier Gas Switching-

- He and H₂ carrier gases were alternated line by line in a sequence table for increased analytical efficiency and higher throughput while demonstrating the flexibility the helium conservation module provides its users while also eliminating manual changes to plumbed gases at the GC EPC Module.

References

1. Agilent GC Calculators and Method Translation Software. Tools available for download from: <https://www.agilent.com/en/support/gas-chromatography/gccalculators>
2. Cristina, M. M.-L. et al. Prediction of Coumarin and Ethyl Vanillin in Pure Vanilla Extracts Using MID-FTIR Spectroscopy and Chemometrics. *Talanta* 2019, 197, 264–269. DOI: 10.1016/j.talanta.2019.01.033.
3. Lingxia, X. et al. Advances in the Vanillin Synthesis and Biotransformation: A Review. *Renewable and Sustainable Energy Reviews* 2024, 189(A). DOI: 10.1016/j.rser.2023.113905